

EM 190 Steady Stand Assist

Product Code: EM-190

The EM190 is a cost-effective and versatile stand aid and transfer support unit, designed to assist in the standing and transferring of users up to 190kg. It encourages more mobile clients to stand un-assisted, making it ideal for those with some mobility but reduced leg strength or balance.

The EM190 is easy to use, featuring legs that open wide to navigate chairs and other obstacles. Its dual braked rear castors keep the unit stable while the patient raises themselves. Once standing, the folding seat pads move into place to prevent falls, allowing for safe, simple transfers between beds, chairs, and toilets without a wheelchair. Suitable for community, aged care, and hospital environments.



KEY FEATURES

High Weight Capacity: Designed to safely assist in the standing and transferring of users weighing up to 190 kg.

Braked Rear Castors: Two rear brake castors provide stability and security as the user raises and lowers themselves from a seated position.

Adjustable Padded Knee Supports: Provides height and tilt adjustment on the padded knee supports for correct, comfortable, and stable user positioning.

Smooth Non-Slip Castors: Large 10 cm non-slip dual rubber castors ensure a smooth, stable transfer and easy manoeuvring on any surface.

Adjustable Opening Legs: Manually opening legs allow the carer to easily navigate wide obstacles such as chair legs or wheelchairs.

Folding Seat Pads: Padded, folding seat pads allow the user to easily step into the aid and can be secured for seated transport.

Secure Hand Grips: Features secure metal handlebars and grab rails, providing safe and stable grip points for both the user and the carer.

Improved User/Carer Interaction: The open design allows for a clear line of sight, promoting better eye contact and communication between the user and carer.

SPECIFICATIONS

- **Length:** Wheels in: 900 mm
- **Length:** Wheels closed: 1000 mm
- **Width:** Wheels in: 600 mm
- **Width:** Wheels closed: 940 mm
- **Turning Cycle:** Wheels in: 900 mm
- **Turning Cycle:** Wheels closed: 1000 mm
- **Product Weight:** 48 kg
- **Max User Weight:** 190 kg / 29 st

